

LATITUDE™ NXT

PATIENT MANAGEMENT SYSTEM

Communicator Setup and Connection Assistance



Call (+61) 1800 528 488 (Australia) or (+64) 0508 200 886 (New Zealand) to get Connected



Get Connected! Call (+61) 1800 528 488 (Australia) and (+64) 0508 200 886 (New Zealand) toll-free for assistance.



Your doctor is prescribing the LATITUDE™ NXT Patient Management System to help your healthcare team collect information about your implanted device and manage your ongoing care.

Call Boston Scientific Patient Service at (+61) 1800 528 488 (Australia) or (+64) 0508 200 886 (New Zealand) for assistance in getting your Communicator set up. Hours are Monday – Friday, 9 a.m. to 5 p.m AEST/AEDT.

Boston Scientific offers comprehensive connectivity options for LATITUDE™ patients. The following connectivity options are available for most LATITUDE Communicators.

Connectivity Option	Hardware	Home Considerations
Ethernet 	LATITUDE Ethernet Adapter	Requires access to high speed Internet modem/router.
Cellular 	LATITUDE Cellular Adapter	Service provided by a network of cellular providers

Please call (+61) 1800 528 488 (Australia) or (+64) 0508 200 886 (New Zealand) and a Boston Scientific Patient Services representative can discuss your connectivity options, provide information about coverage area,*† and assist with the setup of your Communicator.

LATITUDE™ (NXT) Patient Management System

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE.

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device or at www.IFU-BSCI.com.

Prior to use, please review the Directions for Use (DFU) for full operating instructions.

Important Safety Information

LATITUDE™ NXT Patient Management is a remote monitoring system that gives your health care provider access to your implanted device data. The LATITUDE Patient Management system is not intended to assist with medical emergencies.

If you are not feeling well, call your physician or 911. The Communicator does not provide continuous monitoring.

The Communicator is designed to operate on standard telephone lines like those found in most homes. The Communicator may work on other telephone systems, such as Digital Subscriber Line (DSL) and Voice Over IP (VoIP) Internet systems, if those systems provide an analog interface for connecting the Communicator.

The Communicator is designed to work only with the implanted device of the patient for whom it was prescribed. It will not work with other patients' implanted devices and should be used only as authorized by the prescribing physician. The Communicator is not for use with any pulse generator other than a Boston Scientific device. Ask your physician if you have questions about any risks with using the Communicator or your implanted device.

It is very important that the Communicator remain plugged into the power outlet. Your communicator should remain connected to a telephone line, ethernet adaptor or cellular adaptor. Some household appliances and other sources of electromagnetic energy could interfere with the communication between your Communicator and your implanted device. You should be at least 36 inches (3 ft.) away from televisions, VCRs, DVD players, personal computers, and other electronic equipment, when you are using the Communicator.

It is recommended that the customer install a surge arrester in the electrical outlet to which the Communicator is connected. *Rx only.*
92481216 (Rev. A)

Australia and New Zealand: Boston Scientific Pty Ltd | PO Box 332 Botany NSW 1455 Australia
Tel: +61 2 8063 8100 -Fax +61 2 9330 1404

*This program's wireless coverage is limited to the data coverage provided by Vodafone™ and its network of cellular providers. Boston Scientific cannot guarantee coverage.

†Certain terms and conditions apply.

All trademarks are the property of their respective owners.


Advancing science for life™

© 2023 Boston Scientific Corporation or its affiliates. All rights reserved.

CRM-1534107-AA