



# The LUX-Dx™ ICM System\*

## Finding the Right Monitoring Option

The patient app on a mobile device is designed to collect information about irregular heart rhythms stored on your Insertable Cardiac Monitor (ICM) device and automatically sends that information to your healthcare team. You have two options: downloading an app onto your own smart device or using the app on a Boston Scientific provided mobile device.



### Downloadable myLUX™ Patient App ANZ\*\*

Download the myLUX Patient App ANZ on your own smartphone to deliver vital cardiac data to your healthcare team. This option enables you to seamlessly integrate the app into your lifestyle. The downloadable app is free to download and use.

This option may be best if you:

- own a smartphone
- keep your smartphone powered on and near you a considerable portion of time each day
- are comfortable using a smartphone or have a family member who can support you

**If you choose this option, please refer to the back of this page.**



### Boston Scientific provided myLUX mobile device

The Boston Scientific provided mobile device has the myLUX Patient App pre-installed. This device delivers vital cardiac data to your healthcare team.

This option may be best if you:

- do not own a smartphone

**If you choose this option, the mobile device will be provided at your procedure and you will learn more then.**

OR

## Let's get in touch

For additional help please contact the LUX-Dx™ Patient Support  
(Mon - Fri 09:00 to 17:30).

Australia: (+61) 1800 528 488 / New Zealand: (+64) 0508 200 886



## Download the myLUX™ Patient App ANZ before your procedure

### Using your personal phone to get started

1 Open the app store on your mobile device



2 Tap the magnifying glass to search for the "myLUX Patient App ANZ" from Boston Scientific



3 Download the myLUX Patient App



#### Note:

- If you are using an Apple™ device, you may be asked to authenticate with your Apple ID and password, fingerprint, or Face ID to download the app. If you can't remember your Apple ID and/or password, you can go to <https://iforgot.apple.com>
- In the app store if you don't see the **Get button** (on Apple devices) or **Install button** (on Android™ devices), ensure your OS version is up to date and try to download the myLUX app again. If the buttons still aren't visible, your mobile device or OS version is not compatible with this app. Your clinician may recommend the Boston Scientific Mobile Device as an alternative monitoring option.

### Bring the following items to your procedure:

- Your smartphone with the myLUX app downloaded.

If the app is not downloaded, you will also need the following information:

If you have an iPhone™:

☐ Apple ID and password

If you have an Android phone:

☐ Google™ account and password

**Boston  
Scientific**  
Advancing science for life™

\*References to "LUX-Dx ICMs" Includes LUX-Dx™ and LUX-Dx II+™ ICM Systems.

\*\*The downloadable myLUX Patient App ANZ is only available for LUX-Dx II+ICM. It is not available for the LUX-Dx™ ICM.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at [www.IFU-BSCI.com](http://www.IFU-BSCI.com). Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries.

Australia and New Zealand: Boston Scientific Pty Ltd | PO Box 332 Botany NSW 1455 Australia

© 2025 Boston Scientific Corporation  
or its affiliates. All rights reserved.

CRM-2157601-AA