



ImageReady™ MR-Conditional* Model List*

Therapy	Device Name	Device Model Numbers	MRI System Configuration*
ICM ¹	LUX-Dx™	M301	1.5T or 3T
S-ICD ²	EMBLEM™	A209	1.5T
	EMBLEM™ MRI	A219	

The following S-ICD electrodes are labeled as MR-Conditional*

S-ICD ELECTRODES: 3010, 3400, 3401, 3501

Brady Pulse Generators ⁴			
Therapy	Device Name	Device Model Numbers	MRI System Configuration*
Pacemakers	ACCOLADE™ MRI	L310, L311, L331	1.5T or 3T
	FORMIO™ MRI	K289	1.5T
	VITALIO™ MRI	K285, K286, K287	
	INGENIO™ MRI	K185, K186, K187	
	ADVANTIO™ MRI	K085, K086, K087	
CRT-P	VISIONIST™ X4	U228	1.5T or 3T

The following Pacemaker and CRT-P leads and accessories are labeled as MR-Conditional*

Right Atrium and Right Ventricle Leads and Accessories

INGEVITY™ MRI: 7731, 7732, 7735, 7736, 7740, 7741, 7742
INGEVITY+™: 7840, 7841, 7842
Suture Sleeves for INGEVITY™ MRI /INGEVITY™+ leads: 6402
FINELINE™ II EZ: 4469, 4470, 4471, 4472, 4473, 4474
FINELINE™ II: 4456, 4457, 4458, 4459, 4479, 4480
Port Plugs: 7145

Left Ventricle Leads and Accessories

ACUITY™ X4: 4671, 4672, 4674, 4675, 4677, 4678
Suture Sleeves: 4603
Port Plugs: 7148



ImageReady™ MR-Conditional* Model List*

The following ICD and CRT-D leads and accessories are labeled as MR-Conditional*

Tachy Pulse Generators ³			
Therapy	Device Name	Device Model Numbers	MRI System Configuration*
ICDs	RESONATE™	D432, D433	1.5T or 3T
	MOMENTUM™	D120, D121	1.5T
	AUTOGEN™	D044, D045, D046, D047, D174, D175, D176, D177	
	DYNAGEN™	D020, D021, D022, D023, D150, D151, D152, D153	
CRT-Ds	RESONATE™	G447	1.5T or 3T
	MOMENTUM™	G124, G125, G128, G138	1.5T
	AUTOGEN™	G172, G173, G177, G179	
	DYNAGEN™	G150, G151, G156, G158	

Right Atrium and Right Ventricle Leads and Accessories

INGEVITY™ MRI: 7735, 7736, 7740, 7741, 7742
INGEVITY+™: 7840, 7841, 7842
Suture Sleeves for INGEVITY™ MRI /INGEVITY+ leads: 6402
FINELINE™ II: 4479, 4480
FINELINE™ II EZ: 4469, 4470, 4471, 4472, 4473, 4474
Suture Sleeves for Fineline leads: 6220, 6221
Port Plugs: 7145

Left Ventricle Leads and Accessories

ACUITY™ X4: 4671, 4672, 4674, 4675, 4677, 4678
Suture Sleeves: 4603
IS4/DF4 Port Plugs: 7148
ACUITY™ Spiral: 4591, 4592, 4593 (1.5 T only)
EASYTRAK™ 2: 4542, 4543, 4544 (1.5 T only)
IS-1 Port Plugs: 7145 (1.5 T only)

Right Ventricle Leads and Accessories

ENDOTAK RELIANCE 4-SITE™ (DF4): 0262, 0263, 0265, 0266, 0272, 0273, 0275, 0276, 0282, 0283, 0285, 0286, 0292, 0293, 0295, 0296
ENDOTAK RELIANCE™ (DF1): 0170, 0171, 0172, 0173, 0174, 0175, 0176, 0177, 0180, 0181, 0182, 0183, 0184, 0185, 0186, 0187 (1.5 T only)
RELIANCE 4-FRONT™: 0636, 0650, 0651, 0652, 0653, 0654, 0655, 0657, 0658, 0662, 0663, 0665, 0672, 0673, 0675, 0676, 0682, 0683, 0685, 0686, 0692, 0693, 0695, 0696
Suture Sleeves for RELIANCE 4-FRONT™: 6403
DF-1 Lead Port Plug for ENDOTAK RELIANCE™ (IS-1) Leads- Single Coil: 6996 (1.5 T only)

References:

1. MRI TECHNICAL GUIDE IMAGEREADY MR CONDITIONAL ICM SYSTEM 92216913-001 EN Europe 2022-07
2. MRI TECHNICAL GUIDE IMAGEREADY MR CONDITIONAL S-ICD SYSTEM 92346927-001 EN EU 2020-09
3. MRI TECHNICAL GUIDE IMAGEREADY MR CONDITIONAL DEFIBRILLATION SYSTEM 51164993-002 EN AUS 2021-10
4. MRI TECHNICAL GUIDE IMAGEREADY MR CONDITIONAL PACING SYSTEM 97036544-001 EN Australia 2023-04

Disclaimer:

Some products may not be currently approved for commercial supply in Australia and New Zealand, and is for information purposes only.

*When conditions of use are met.

Visit <http://www.bostonscientific.com/imageready> for additional information cardiology/radiology checklists, patient resources, and the MRI technical manual.



2024 Copyright © Boston Scientific Corporation or its affiliates. All rights reserved.

CRM-1889601-AA

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device or at www.IFU-BSCI.com. Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. This material not intended for use in France.
Boston Scientific Pty Ltd | PO Box 332 Botany NSW 1455 Australia | Tel +61 2 8063 8100